

CYYR1-cysteine/tyrosine-rich 1 Gene

Reference: PTXBC036761

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf95; cysteine and tyrosine-rich protein 1; cysteine and tyrosine-rich protein 1 isoform 1,2,2bis,3,4; cysteine and tyrosine-rich protein 1 isoform 1,2,3,4b; cysteine and tyrosine-rich protein 1 isoform 1,2,3b; cysteine and tyrosine-rich protein 1 iso		

For research use only.

GENE DESCRIPTION	cysteine/tyrosine-rich 1
GENE ID	116159
ORIGIN SPECIES	Human
ACCESSIONS	BC036761
VECTOR	pDONR223
NCBI SYMBOL	CYYR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgctccgaggtaccgtgctccagggtcttctccgaagtgtgctcttctgtacgcagatgattgcttgcctcagtggtgcaagattgcaaatcttactgctgtgatggaaccacgcctactgttgcctactacgcttatattgggaatatctctcgggcactgcaat
tgcgggcattgttttggaaatagttatcatggtggtcattgctgggattgccatgcatctgcatgcatgaagaaccacaggcgacccgctggtggcatctcaggacgactcacatcaacacgctctctctatctggaaccacccctacggtcagcaccagatggaaat
tgtgcagactgctctcctactccccacccacagggtccagcacagcgttctccacccctcttctcctggaacgcaaggaaataa