



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NEGR1-neuronal growth regulator 1 Gene

Reference: PTXBC036771

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

DMML2433; IGLON4; KILON; Ntra; neuronal growth regulator 1; IgLON family member 4; a kindred of IgLON; neurotractin

For research use only.

GENE DESCRIPTION	neuronal growth regulator 1
GENE ID	257194
ORIGIN SPECIES	Human
ACCESSIONS	BC036771
VECTOR	pDONR223
NCBI SYMBOL	NEGR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggctgcactctaactgtgcgaagcttcctcctaagatfatatgacatctcaaatgatgacgcgcfaatgaaggaaccaacgcctcttactctgttggccactgggaacacagagccttcattcttggcgacacatctccccacgaaaaacatttggaaatggacaatatttggacattat
ggaaatcaacaagggacacaggctgggggaatatgaatgcagtcgggaaatgatgtgtcattcccgatgtggaggaaatgaagaattgtgtgacacttggcctctacttactcaggaaataaatctggcaccgtgacccccgacgcagtgccgtgataagatgtgaaggtgaggtgtgtgccgcct
cacgccttgaatgtgtacaaagagagagaagcgtcttcaattgtgtcccaacgaagaattattatcaaaattttagcaacagatcaacttcaactgttcaacgcgtgacacaggaagcatttggcgaaattactctgtgtgcttgcacacagtagtgccacaacaaatgcgagctgcctcttaacc
tccaaatgcacgcgaatctgaattgaattgagagagagcctcttgcctctctctctctctctctgtgttgaactcactctctctcttcaacttctactctctcaagaattgccattctcaataa