

CRYBB1-crystallin, beta B1 Gene

Reference: PTXBC036790

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CATCN3; CTRCT17; beta-crystallin B1; beta-B1 crystallin; eye lens structural protein; crystallin beta B1		

For research use only.

GENE DESCRIPTION	crystallin, beta B1
GENE ID	1414
ORIGIN SPECIES	Human
ACCESSIONS	BC036790
VECTOR	pDONR223
NCBI SYMBOL	CRYBB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcagctgcaaggcctcgccctggccacagtggtgaaccagggcctgacaccaaggggaagggggcccccacctgcaggaaacatccctagtcggcactaccctggcccaacaacgtgcctattaccagcgccaaggcgagggaactgctctctgggaact
acaggctgggtgtcttgaactggaaactccaggccgtcgagcagaattctcggggagtgctcaaatctggcagaccgtggtcttcgaccgtgtgcgcagcatcattgtctcggcgggaccctgggtgcctttgagcagtcgaacttcgcggggagatgtcatctggagaagg
gcgagtagccctcgtggaacacatggtcgagcagctaccgagtgatggctcatgtctctccggccatcaaatggatgccaggagcacaataatctccctgtttgaagggggccaacttcaagggcaacacatagagatccagggggacgacgacccagctctctgggtctacggc
ttcagtgaccgctggggcagcgtgaaggtctccagtggaacatgggttggctatcagatctggctaccgcccgtaccagtagctcctagagcctggtgacttcggcactgggaatgagtgaggagccttcagccacagatgcagtcctcgtcgtcgtgacaagcagtgccac
ctcgaggggtccttcctgtcctggccacagagcccccaagtga