

C1orf131-chromosome 1 open reading frame 131 Gene

Reference: PTXBC036800

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf131; chromosome 1 open reading frame 131		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 131
GENE ID	128061
ORIGIN SPECIES	Human
ACCESSIONS	BC036800
VECTOR	pDONR223
NCBI SYMBOL	C1orf131
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcaggagcaagggccgggtctccacgctccagttctccgacactcttgacgctctgctccagaacctttacgactttgaggtacagaagtgaaacagacagaagaagatcataaagaaaaggaaaacaagaagagagatgtgatggcttcagcggttcggcagc
agagccatctccctaccctgggtctctcataagggccagaggaagagcgcttcgagctcttcaaggaacttagagaagagcggttcgctctctggtgacccccacaggaccagagatccttgcctgcgagttcctcctctccctaaagaacaatagggaacaagttagaagtgtg
agaatttcacagcaataaaaaaagaaattgacgccagatcataaacaagaacacaaaggcctaatacctagttgtttggagagagatgtggatacacaagaatttaacctagaaaaagctcgtttagaagtgcaccggtttggtatcacgggttatgaaaaggaaaaggagagaatcctggac
agggaacgtgccattatgctggcgctaaagcctcctaaaaagagttatgtgaattacaaggtttttacaggagcaaataaagaaaaaaggcagcaagggaagaaaaagagactggcccaagaacagatatttcaagaaaaagaaggaaggacaggagacaggaaatccaa
aaagaagtcgctccagatatttgcataatggacgattggacaggttgaaaaattcaaaaatggaacactgattctgagccagttgatatcaagaaaaataattcttcagagtggtcccaatga