



Reference: PTXBC036812

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GALNTL1; GALNACT16; GalNAc-T16; polypeptide N-acetylgalactosaminyltransferase 16; UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase-like protein 1; UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase 16; UDP-N-acetyl-al		

For research use only.

GENE DESCRIPTION	UDP-N-acetyl-alpha-D-galactosamine:polypeptide acetylgalactosaminyltransferase-like 1	N-
GENE ID	57452	
ORIGIN SPECIES	Human	
ACCESSIONS	BC036812	
VECTOR	pDONR223	
NCBI SYMBOL	GALNTL1	
LENGTH	ProteoGenix bp	
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG	

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr