



Reference: PTXBC036827

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD85D; ILT-4; ILT4; LIR-2; LIR2; MIR-10; MIR10; leukocyte immunoglobulin-like receptor subfamily B member 2; CD85 antigen-like family member D; Ig-like transcript 4; leucocyte Ig-like receptor B2; leukocyte immunoglobulin-like receptor, subfamily B (with		

For research use only.

GENE DESCRIPTION	leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 2
GENE ID	10288
ORIGIN SPECIES	Human
ACCESSIONS	BC036827
VECTOR	pDONR223
NCBI SYMBOL	LILRB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]