

GTPBP8-GTP-binding protein 8 (putative) Gene

Reference: PTXBC037163

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC135; GTP-binding protein 8; GTP binding protein 8 (putative)		
For research use only.			

GENE DESCRIPTION	GTP-binding protein 8 (putative)
GENE ID	29083
ORIGIN SPECIES	Human
ACCESSIONS	BC037163
VECTOR	pDONR223
NCBI SYMBOL	GTPBP8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccgggctgctgggctggagcggggaagactcttgaatgcctgcgtgctagagcagctgagccgtataatgacgtcccaagcttctgaggtgctgctggcgaagcagcagctgaggaagctgtacccgtgcaggagtagagcgggttcctgccc
cctacgggaggcgaagacctcactgcgtatcttgaaccaagccggaggacatagccaggcgggacaacatcttcacggccactgaacggaacccgacgactacgtcagctccgctccgctatcgaccacgccccggaccttcgcgccagagggtgtttatagggcagaag
caatgttgaaaatcatcttaatacaaggctttatctactggccctgaggttgagtcagagctctccaaaaaccaggacacacaaagaaatgaatttttcaaaagttggaaaacattttacagtggtggacatgccaggttatggctttagagccactgaagattttgtgacatgtagaga
cctatctaaaagaacgaaggaaactgaagagaacattttatagtgatagcgtgttggaattcaaaaaacagacaatttgcctatgaaatgtgtgaagaatttgacattacattatgtgattgtattaacaaaattgacaattctccaaggacatctttaaaacaagtgcttcagatccagaa
attgttaacatgaaaactcaaggatgtttctcagttgttctcgtgaagtgtgtgacctttctggaatccacctgttgagatgctttatagccaggtgaacaggaggtcttgactaa