



Reference: PTXBC037283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FACE-1; FACE1; HGPS; PRO1; STE24; Ste24p; CAAX prenyl protease 1 homolog; farnesylated proteins-converting enzyme 1; prenyl protein-specific endoprotease 1; zinc metallopeptidase STE24 homolog; zinc metalloproteinase Ste24 homolog; zinc metallopeptidase S		

GENE DESCRIPTION	zinc metallopeptidase (STE24 homolog, <i>S. cerevisiae</i>)
GENE ID	10269
ORIGIN SPECIES	Human
ACCESSIONS	BC037283
VECTOR	pDONR223
NCBI SYMBOL	ZMPSTE24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]