



Reference: PTXBC037308

|                |                                                                                                                                                                                             |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | YY1 transcription factor; transcriptional repressor protein YY1; DELTA; INO80S; NF-E1; UCRBP; YIN-YANG-1; INO80 complex subunit S; YY-1; Yin and Yang 1 protein; delta transcription factor |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | YY1 transcription factor                                                                          |
| <b>GENE ID</b>          | 7528                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC037308                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | YY1                                                                                               |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]