



Reference: PTXBC037556

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA0182; CRHSP24; genetic suppressor element 1; CTC-786C10.1; Gse1 coiled-coil protein homolog; Gse1 coiled-coil protein		

For research use only.

GENE DESCRIPTION	KIAA0182
GENE ID	23199
ORIGIN SPECIES	Human
ACCESSIONS	BC037556
VECTOR	pDONR223
NCBI SYMBOL	KIAA0182
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]