

TPM4-tropomyosin 4 Gene

Reference: PTXBC037576

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-108; tropomyosin alpha-4 chain; TM30p1; epididymis secretory protein Li 108; tropomyosin 4		

For research use only.

GENE DESCRIPTION	tropomyosin 4
GENE ID	7171
ORIGIN SPECIES	Human
ACCESSIONS	BC037576
VECTOR	pDONR223
NCBI SYMBOL	TPM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggcctcaactccctggagcgggtgaacgcaagatccaggccctgcagcagcggcgacgagggcggaagaccgcgcgagggcctgcagcgggagctggacggcgagcgcgagcggcgagaaagctgaaggtgatgtggccgcccctcaaccgacgcacccag
ctcgttgaggaggaggttgacagggcctcagggaacgactggccacggccctgcagaagctggaggaggcagaaaaagctgcagatgagagtgagagaggaatgaaggtgatagaaaacggccatgaaggtgaggagaagatggagattcaggagatgcagctcaagaggc
caagcacattgcggaagaggctgaccgcaaatcaggaggtagctcgttaagctggtcatctggagggtgagctggagaggcgagaggagcgtgcggagggtgctgaactaaatgtggtgacctggaagaagaactcaagaatgttaactaactctgaaatctctggaggctgca
tctgaaaagtattctgaaaaggaggacaatatgaagaagaattaaacttctgtctgacaaactgaagaggctgagaccgtgctgaatttcagagagaaacgggtgcaaaactggaagagacaattgatgacctggaagagaaactgccaggccaaagaagaacgtgggctta
catcagacactggatcagacactaaacgaacttaactgtatataa