



Reference: PTXBC037776

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NTR2; neurotensin receptor type 2; NT-R-2; levocabastine-sensitive neurotensin receptor; neurotensin receptor, type 2; neurotensin receptor 2		

For research use only.

GENE DESCRIPTION	neurotensin receptor 2
GENE ID	23620
ORIGIN SPECIES	Human
ACCESSIONS	BC037776
VECTOR	pDONR223
NCBI SYMBOL	NTSR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1