

LOC153364-similar to metallo-beta-lactamase superfamily protein Gene

Reference: PTXBC038230

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to metallo-beta-lactamase superfamily protein
GENE ID	153364
ORIGIN SPECIES	Human
ACCESSIONS	BC038230
VECTOR	pDONR223
NCBI SYMBOL	LOC153364
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgctcagtggtacgccacaagtctctagggcatggtatcttctggaatcaagaacgtttctacagtcgggcaaccgtgccacatctggctggcgccgctccgagcaggacgtgggtgatcgatacaggcctggggctgcgcagcctccggagtaacctgtactcctcg
gcctcttgaggaccgagagggccaaagagagcgcgcgcggccactgcttgcctggccaccacgtgcacttcgaccactcggcgccctaccagttcgaccgctggcagtgaccacgccgagggcgctgctcgcgggacaactttgagaccgtgacctg
gcttccgatacgaggtgggtggcgccagccccggctggaggccagacagttccgggtacagggcggtgcagccaccctcatctcgaggtgggatggtgacacacttgggtgacagacagctcactgttatgcacatgccaggtcactccagggcgatgttcttcatg
acaaagaccgaaagattctctcagtgagacgtcgtgtatggtgacactgattgactggctcccatacagcaggataagtgactatgttggaaactgtgaacgtctaatagaattagtgacagaggtctggtagagaaggtgcttctgggcacttcaataccttgggtgctgaaagcctt
ttcgattggcttcaactatattcaaaagctgggataatgcacaaagtcttacttttgcctatgcgatctcttgaagtttagctctacgtgtaacaaattctaggacctgcctag