



Reference: PTXBC038237

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SLEB6; integrin alpha-X; CD11 antigen-like family member C; complement component 3 receptor 4 subunit; integrin alpha X; integrin, alpha X (antigen CD11C (p150), alpha polypeptide); integrin, alpha X (complement component 3 receptor 4 subunit); leu M5. al		

GENE DESCRIPTION	integrin, alpha X (complement component 3 receptor 4 subunit)
GENE ID	3687
ORIGIN SPECIES	Human
ACCESSIONS	BC038237
VECTOR	pDONR223
NCBI SYMBOL	ITGAX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]