



Reference: PTXBC038239

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JTK14; TIE; tyrosine-protein kinase receptor Tie-1; tyrosine kinase with immunoglobulin and epidermal growth factor homology domains 1; tyrosine kinase with immunoglobulin like and EGF like domains 1		

For research use only.

GENE DESCRIPTION	tyrosine kinase with immunoglobulin-like and EGF-like domains 1
GENE ID	7075
ORIGIN SPECIES	Human
ACCESSIONS	BC038239
VECTOR	pDONR223
NCBI SYMBOL	TIE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

More info: Please contact us for any question: contact@proteogenix.fr