

ELOVL4-elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 4 Gene

Reference: PTXBC038506

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3-keto acyl-CoA synthase ELOVL4; ADMD; CT118; ISQMR; SCA34; STGD2; STGD3; elongation of very long chain fatty acids protein 4; ELOVL FA elongase 4; cancer/testis antigen 118; elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 4;		

For research use only.

GENE DESCRIPTION	elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 4
GENE ID	6785
ORIGIN SPECIES	Human
ACCESSIONS	BC038506
VECTOR	pDONR223
NCBI SYMBOL	ELOVL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggctcctggactcggagccgggtagtgcctaaacgtagtgtccacggcactcaacgacacggttagagttctaccgctggacctggccatcgcagataagcgtgtggaaaattggcctctgatcgcggtctcctggcctacactaagataagcactctttatctcctgtttgtgtgct
gggtccaaaatgatgaaggaccgagaaccltttcagatgcgtctagtctcattatctataattttgggatgttttgccttaacctctttatctcagagagttatcatgggatcatataatcgggatatactatattgccagagtgtggattattctaataatgttcagagtcaggatagctgct
gctctgtgtgtgactttgtatctaaaggagttgagtttggacacagtggtttttattctgagaagaaaaacaacaaagtttcttctcctatgtatcatcactgtacgatgtttacctgtgtgtgattggaattaaagtggtgtcaggaggacaaagcattttttggagccagttgaattccttat
ccatgtgattatgtactacatactgggttaactgcatttggccatggattcagaataatcttgggtgaaacgatacctgactatgttgcaactgatcaattccatgtgaccattgggcacacagcactgtctctttacactgactgcccttcccaaatggatgcactgggcttaattgcctatg
caatcagcttcattttcttcttaactctacattcggacatacaaaagagcctaagaacaaaagctggaaaacagccatgaatgtgttgcagcaaatggtgtgagcaaatcagaaaaaacatcatgataaaaatgaaaaagcagaaaatgaaaaagcaaaaggagattaa