

C12orf60-chromosome 12 open reading frame 60 Gene

Reference: PTXBC038836

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C12orf60; chromosome 12 open reading frame 60		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 60
GENE ID	144608
ORIGIN SPECIES	Human
ACCESSIONS	BC038836
VECTOR	pDONR223
NCBI SYMBOL	C12orf60
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcagagtcagaaaaggataaaagagagactgattcaagctgccaatgtctcttcacgtacaaagatctgtctctgtcataaacacactgactgaattgtttagccgcagatgaataactcaaatctttgatggctgtgagaaacaatagttacattaaggattttttgagcaaatgctc
aaaatttttaaggagatgcaatctgtatgtgcaagacatgacaaaatcaaaaggagctctttatgttccaaggttgcaatggccatgtgctctgtgttcagaagagtaccaatgtagaggagttgcatcagtcagctaaagaagtattcaaaagtcccatacggcagtcacatctctgtg
ctaaacagcagtaacatccttgggagctggaatctctcttcacacttgatgaattcccacatgaatcttcaatttaagtgacttctatacagaagacaccaagagcaatcagatgtcaccacatctgagagaaccagaagtcctccaggttctccaaaaccactatgatacacacttga
aaaaactgcaggtgtactaaaaactgaggattccaaaatcccacaaagtcagcagcagattgttggacaattgtcaaggctatgggaccaatcttagagatctccaaaaagcgataaagactatggaatgaatatttctgtgtttaagaaagccagtgacaagtag