

GJB6-gap junction protein, beta 6, 30kDa Gene

Reference: PTXBC038934

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CX30; DFNA3; DFNA3B; DFNB1B; ECTD2; ED2; EDH; HED; HED2; gap junction beta-6 protein; connexin 30; ectodermal dysplasia 2, hidrotic (Clouston syndrome); gap junction protein, beta 6, 30kDa; gap junction protein beta 6		

For research use only.

GENE DESCRIPTION	gap junction protein, beta 6, 30kDa
GENE ID	10804
ORIGIN SPECIES	Human
ACCESSIONS	BC038934
VECTOR	pDONR223
NCBI SYMBOL	GJB6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattggggacgctgcacacttcatcgggggtgtcaacaacactccaccagcatcgggaaggtgtgatcacagtcatttatttccgagtcgatcctcgtgggtgcccaggagtggtgggtgacgagcaagaggacttctgtgcaacacactgcaaccgggatgca
aaaatgtgtctatgaccacttttcccggtgtccacatcggcgtgtgggcccctcagctgatcttctcaccaccagcgtcgtgggtgcccagtcgtgtggtgacctacagggcacgaaaccactcgcaagttcaggcgaggagagaagagggaatgatttcaagacatagaggacatt
aaaagcagaaggttcggatagaggggtcgtgtgtgtgacgtacaccagcagcatctttccgaatcatcttgaagcagcctttagtattgtgttttacctcttacaatgggtaccacgtgccctgggtgttgaatgtgggtgacccctgcccccaacctgttgactgcttttctagg
ccaacagagaagaccgtgtttaccattttatgttctgctgtgtgttgcgtgcttaacgtggcagagttgtgctacgtgctgaaagtgtgttttagagatcaaaagagacagacgcaaaaaatcacccaatcatgccctaaaggagagtaagcagaatgaatgaatga
gctgatttcagatagtggtcaaaatgcaatcacaggtttcccaagctaa