

CLEC1A-C-type lectin domain family 1, member A Gene

Reference: PTXBC039072

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLEC-1; CLEC1; C-type lectin domain family 1 member A; C-type lectin-like receptor-1		

For research use only.

GENE DESCRIPTION	C-type lectin domain family 1, member A
GENE ID	51267
ORIGIN SPECIES	Human
ACCESSIONS	BC039072
VECTOR	pDONR223
NCBI SYMBOL	CLEC1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggccaagtacagcagcagcaggacatgctggatgatgctgggacaccacatgagcctgattctcaaggctctgccacaactcggcatccagagccccggcgacagagcagcagggtccctcttcaagtgccgaccagtggccctgacctgctgactttgtgcttgggtgctgctgatagggctggcagccctgggcttttttttctagctaccagctctccaatactgggtcaagacacatttctcaaatggagaagattaggaaatagctcccaagagtgcaatctctcaagtcagaataaagctgcagggaagtctgcagcatgtggctgaaaaactctgtcgtgagctgtatacaaaagctggagcacacaggtgcagccctgtacagaacaatggaatggcatggagacaattgctacagttctataaagacagcaaaagtgggaggactgtaaatatttgccttagtgaaaactctaccatgctgaagataaacaagaagacctggaattggccgctctcagagctactctgagtttttctactcttattggacagggcttttgcgccctgacagtggaaggcctggctgtgatggatggaaccccttctactctgaactgttccatattataatagatgtcaccagccaagaagcagagactgtgtggccatccttaattggatgatctctcaaggactgcaagaattgaagcgttgtgtctgtgagagaaggcaggaatgtggaagccagagagcctccatgttccccctgaacattaggcgaagggtgactga