



Reference: PTXBC039579

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OBPH1; ORP5; oxysterol-binding protein-related protein 5; ORP-5; OSBP-related protein 5; oxysterol-binding protein homolog 1; oxysterol-binding protein homologue 1; oxysterol binding protein like 5		

For research use only.

GENE DESCRIPTION	oxysterol binding protein-like 5
GENE ID	114879
ORIGIN SPECIES	Human
ACCESSIONS	BC039579
VECTOR	pDONR223
NCBI SYMBOL	OSBPL5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr