

C11orf53-chromosome 11 open reading frame 53 Gene

Reference: PTXBC039669

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C11orf53; chromosome 11 open reading frame 53		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 53
GENE ID	341032
ORIGIN SPECIES	Human
ACCESSIONS	BC039669
VECTOR	pDONR223
NCBI SYMBOL	C11orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagggtcgcgggtcacgtcaggttactacgggtgcagaagatcttcttctgactcagactccacacagtaaacagttttcaaatgacgtctacacccagcgtggggaagccgtttccctgtgagtcctccgagggcagagccatgcggctctctggagccctacttcccca
ggagccctacggagactaccggcctccggcgtgacgccaacggggctctctgttcagcgctcgcctaccgcccctcctgccaccgcccctcccgagaccagctcacttctttagggactcatgggagcagacgttgcctgacggtctcagccagcctgacctgtgtc
tgccgatccctgtgaccttgccaccagcagagttgcctctccagcttgagtcgggagcatcggccagcacagggctcaagctgggggtcaccctggctggggtcagtcatactcgtgcatgctctggaagatctgcaccaactccgggtaccctaccctccctcttac
ccctcacccttcatgacggtgtcaaatgacctaccgccaaggtggggcactctccagatgaggagcagacaccgggtccctccatgaccttcccttgggtgaaagaagatgggagttatgcctgggggtcatatgaatgccgagagcttattga