

TRDV2-T cell receptor delta variable 2 Gene

Reference: PTXBC039714

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	T cell receptor delta variable 2
GENE ID	28517
ORIGIN SPECIES	Human
ACCESSIONS	BC039714
VECTOR	pDONR223
NCBI SYMBOL	TRDV2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagagatcctccctcatctctctctctctggcaggagtcgagccattgagttggcctgaacacaaacagtcgctgtgcaatagggtccctgccacccctcagggtgctccatgaaaggagaagcgtcggtaactactatcaactgggtacaggaagacccaa
ggtaacacaatcatttcatataccgagaaaaggacatctatggccctgggttcaagacaattccaaggtgacattgataatgcaaaagacctggctgtacttaagatactgcaccatcagagagagatgaagggtcttactactgtgctgtgacggcagaggactggggatacagatt
caggaaacacacacctctgtcttggaaaggcacaagacttctgtgattgcaaatatccagaacctgacctgcccgtgtaccagctgagagactctaaatccagtgacaagctgtgtcctattcaccgatttgcattctcaacaaatgtgcacaaagtaaggattctgattgtatatacaca
gacaaaactgtctagacatgaggtctatggacttcaagagcaacagtcgtgtggcctggagcaacaaatctgacttgcattgtgcaaacgccttcaacaacagcattattccagaagacaccttctcccagccagaaagtctgtgtgctcaagctgggtcgagaaaagctttgaaaca
gatacgaacctaaacttcaaaacctgacgtgattgggtccgaatcctctcctgaaagtggccgggttaattctgctcatgacgtcggtgtgtgtccagctga