

TMEM106B-transmembrane protein 106B Gene

Reference: PTXBC039741

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 106B		

For research use only.

GENE DESCRIPTION	transmembrane protein 106B
GENE ID	54664
ORIGIN SPECIES	Human
ACCESSIONS	BC039741
VECTOR	pDONR223
NCBI SYMBOL	TMEM106B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaagtcttctcatttgcctttgcaagcaagaagatgcttatgatggagtcacatctgaaacatgaggaatggactggtaatagtgaagtcataatgaagatggagaatggagatgtctcagttccatatgtggaatttacaggaagagatagtcacctgcccta
ctgtcagggaacaggaagaattcctagggggcaagaaaacaaactgggtgcatgattccatagtgatcagagattaagggcaagaagaacaaagctgtatgtggtctctgtgtttgtctgtctactcttctggattggctgtgttttcccttccctcgtctatcgacgtgaaatacat
tggtgtaaaatcagcctatgtcagttatgtatgagcgtacaatttatttaaatcacaaacacactaaataatacaaaacaataactattactctgtcgaagttgaaaacatcactgcccaagttcaattttcaaaaacagttattggaaggcagccttaaacataagcattattggtcca
cttgatatgaacaaattgattacacagtacacctgtatagcagaggaatgagttatgtatgatttctgtactctgatatccaaagtgcaatacatgactcatgcaagttactgtgacaacaacatacttggccactctgaacagatcccaggagaggtatcagtatgtcgac
tgtggaagaacacaaactatcagttgggcagtcgtgaatttaaatgtacttcagccacaacagtaa