



Reference: PTXBC040007

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	prosaposin receptor GPR37; EDNRBL; PAELR; hET(B)R-LP; ETBR-LP-1; G protein-coupled receptor 37 (endothelin receptor type B-like); Parkin-associated endothelin receptor-like receptor; endothelin B receptor-like protein 1; G protein-coupled receptor 37		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 37 (endothelin receptor type B-like)
GENE ID	2861
ORIGIN SPECIES	Human
ACCESSIONS	BC040007
VECTOR	pDONR223
NCBI SYMBOL	GPR37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr