

IQCF2-IQ motif containing F2 Gene

Reference: PTXBC040047

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IQ domain-containing protein F2; IQ motif containing F2		
For research use only.			

GENE DESCRIPTION	IQ motif containing F2
GENE ID	389123
ORIGIN SPECIES	Human
ACCESSIONS	BC040047
VECTOR	pDONR223
NCBI SYMBOL	IQCF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggggttcgattttgtaccaaaggcaattttaattttgtttataattgaggtatgttgaagaaagcattgaatggaagacattgcagaagaagaacagcagaaaaactagaataagaacaaagcagctgtaagatccaggcctgggtggcggggcacctgggtgcgcag
gacactgctgcatgcagccctcaggcctggataattcagtcgtgtggcggatgacgctgtcaggggtgctggagaagaacggcaggcagctctgacgctacgcaaccagagagagggcagtgatcaagctccagtccttgggtccgtatgtggcgtgtccgtggcgatactgcc
agggtgctcaatgccatctacatcatccaggccactggcaatgccacaactgccagactgcgctctcctccaggccactgtgtgtgtcacagccactacctgcagttccacattgagatcatcaactcctaa