

SNHG11-small nucleolar RNA host gene 11 (non-protein coding) Gene

Reference: PTXBC040237

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A930034L06Rik; AI854265; AL022637; AW319638; E130013N09Rik; small nucleolar RNA host gene 11 (non-protein coding); small nucleolar RNA host gene 11		

For research use only.

GENE DESCRIPTION	small nucleolar RNA host gene 11 (non-protein coding)
GENE ID	128439
ORIGIN SPECIES	Human
ACCESSIONS	BC040237
VECTOR	pDONR223
NCBI SYMBOL	SNHG11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggggcaggcccgaggagtggtcttccaagaacctgggtggcctccaaggccggtgtgtgtacctcctccacgacaaaagggaactgaggcccgaggggagtgaggagagccggctggacgtcagggccagccgctggtgcagtggtcgtccctctgccgg
gggtggccctcggtttcgcgtgtcctcgggaaagagactggcgggtgagccgcgcctcggccttcgtgggctaa