



Reference: PTXBC040291

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM116A; FAM116A; AFI1A; protein DENND6A; DENN domain-containing protein 6A; DENN/MADD domain containing 6A; family with sequence similarity 116, member A; DENN domain containing 6A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 116, member A
GENE ID	201627
ORIGIN SPECIES	Human
ACCESSIONS	BC040291
VECTOR	pDONR223
NCBI SYMBOL	FAM116A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr