

MCFD2-multiple coagulation factor deficiency 2 Gene

Reference: PTXBC040357

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	F5F8D; F5F8D2; LMAN1IP; SDNSF; multiple coagulation factor deficiency protein 2; neural stem cell-derived neuronal survival protein; multiple coagulation factor deficiency 2		

For research use only.

GENE DESCRIPTION	multiple coagulation factor deficiency 2
GENE ID	90411
ORIGIN SPECIES	Human
ACCESSIONS	BC040357
VECTOR	pDONR223
NCBI SYMBOL	MCFD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccatgagatccctgtctcagaaccccttctgtgtgctgtctgtggcctttgtgccccaggcgccagggtgaggagcctgcagccagcttctcccaacccggcagcatgggcctggataagaacacagtgacgaccaagagcatatcatggagcatctagaaggtgtcatc
aacaacccagaggcggagatgtcgcacaagaattgcagctccattctcaaatgcatgattatgatggcaataattgcttgatggcttagaactctccagccatcactcatgtccataaggaggaaggagtgaaacaggccaccactaatgagtgaagatgaactgattaacataata
gatggtgtttgagagatgatgacaagaacaatgatggatcatgtactatgctgaatttgcaaaatcactgcagtag