

## CEP120-centrosomal protein 120kDa Gene

Reference: PTXBC040527

|                |                                                                                                                                                 |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CCDC100; SRTD13; centrosomal protein of 120 kDa; centrosomal protein 120kDa; coiled-coil domain-containing protein 100; centrosomal protein 120 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | centrosomal protein 120kDa                                                                        |
| <b>GENE ID</b>          | 153241                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC040527                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CEP120                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggtctccaaatccgacaaattgctcatcgtcgtgtccatcctagaaggtcgcatttcccaaacgtccaaagcatatgctttagtggagcaagtttgatggagaacagttggctactgacctgtggaccacactgaccagccagaattgctactgagttgcttgggaaattgacag  
gaaagcgttcatcagcacagggtacagcgtactcctatcaaaactccaatgttttgccttggaactcgttaacttcagccaaggaaacataggttacatcgttctggatttaagaaccgtcagaacaacaaagcaggcaccaaaatggaccagttgctgagtaataatacaccaaatcaagt  
ctgagatacagataaagtattgcttggaaaccgatacaaaagccaccagtggtatagctttaaagcaaaaggggctccccctcgagatggaaaagtaacctgccaatcctggctggacttgacccaggacattggctgtgctgaatgaaggagggtaccatcagattggaccagcag  
aatactgtactgactcctttatgtcagtgaccatagcatttgcacccagttggaacagtttaattccatgtaccatgaaactccagaagacagcctgagttttctttactactctttactgggaaatgatgttacaatgaaccttcaatgatttaacacccaactttgagccagagagag  
catcagttcgatcctgtagcagttgtagaaattctcgtttacctggctctcagttcaactgcagattcacctctgctgtggagaccagtcacttggaaagacagaataacctttaactggattacttaaaaaggcagtacagaatcaaccagcaccagtcacagtcgaaggtcctttac  
ccttgaccctccaacagagccaacagaagcttgacatttctgtggagctagcccaactgtggaggtgtctgtggctctgcagagagaagcatagactccaggatgcatttggattcagctctggatcatatttctctgttctctctttctgtcctagatgctataaggaa  
atttgcataattacgaagaaatag