

C22orf30-chromosome 22 open reading frame 30 Gene

Reference: PTXBC040859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf30; protein PRR14L; proline rich 14-like protein; proline rich 14 like		
For research use only.			

GENE DESCRIPTION	chromosome 22 open reading frame 30
GENE ID	253143
ORIGIN SPECIES	Human
ACCESSIONS	BC040859
VECTOR	pDONR223
NCBI SYMBOL	C22orf30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagagctcttcatgacccaggtttacacagggcctgaaaggggttacggctccagcctccatagcagacaaggcttctgtctctgctcctactcgggtggcagagtcctatccattggagccagcatggctctctgtctctccttgcgaatctcttcttcttccctcactgcaagcg
gcaaccaagctctgggaccacaagcagccacacattgttaccatattgtcctcttccagcctggaagctacatataacacagcggcagtcagacgagcctggagcctccattccctgcttgggtacaaagcttctgttggtagcagaatcagctgtcagcaagcctctgttcaagcctc
tgagttccaggttcgtgattgatgagctggatgggtgaaagcagcatgcccctgccacagagcagcccccaacagaagagggctgagccagagaagaggccaaagaaagtctcacagattcgcattccgaaaaccattcctaggccagatcctaatttaccctcatgggc
cttctcgacccaaaaggttaagaagaaggagtttagtttagaagagatataacacaagaattataatctctcctcgaacaggggacgtggattctagctgcccagagcactgacagcgaaacgggcacttctgttattcggatattgcgagcgccaagcgagccgacccact
gctggcctcacgcgtctgattggctgcagctctccggaggggcggggccatccgagaagctatccatttgaggaggctgctagtccttctgtgttacctgccggctggagagcggggagggtcacgtgtgtaagactgccacgtgggcccggcggtgactcaagcctgtaa