

SYT16-synaptotagmin XVI Gene

Reference: PTXBC040924

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHR14SYT; SYT14L; Strep14; syt14r; yt14r; synaptotagmin-16; CTD-2277K2.1; chr14 synaptotagmin; synaptotagmin 14-like protein; synaptotagmin XIV-related protein; synaptotagmin XVI; synaptotagmin 16		

For research use only.

GENE DESCRIPTION	synaptotagmin XVI
GENE ID	83851
ORIGIN SPECIES	Human
ACCESSIONS	BC040924
VECTOR	pDONR223
NCBI SYMBOL	SYT16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccgagagagaatgatgggagagaaattctatctcagccacctgcaccagaaggggaaatgaaagtgactctggtctggagccaagaagtaataataagcagtgagggtctccgctcagcccatctcggtttctcacagtatagttacttccacgcagtcgtgtctcat
ggagggggccagagctgttgggtgggctctgtacaatgccacaacgggctgattatctgtggaatgatcaaaaggcagccatttccgaaacctgctgttaaccgagcacctgatacatatgaaaaacttttctcctaattctgtgggtcaagagatgtcccgttgcaagacgtccattc
ggcgtggtcagcccaatcctgtctataaggagaccttgggttccaggtggccctctttcagctgtctgatgtcacgttgatgattccggtttataacaggcgtactatgaagcgtaaagagatgattggctggatggccctgggcccagaacagcagtgaggaggagaacaagatcactggga
ggagatgaaggaaaccaaggccagcagatctgcagatggcacactttgctggaatcctag