

TCL6-T-cell leukemia/lymphoma 6 Gene

Reference: PTXBC041075

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TNG1; TNG2; T-cell leukemia/lymphoma 6 (non-protein coding)		

For research use only.

GENE DESCRIPTION	T-cell leukemia/lymphoma 6
GENE ID	27004
ORIGIN SPECIES	Human
ACCESSIONS	BC041075
VECTOR	pDONR223
NCBI SYMBOL	TCL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagccagagtgacacagaggagggccctggatggatgcattgggaagatcacgggaattacaagtacatcctgaagtatgacacaaatgttcaaatgaagcctctgctaagttccagaggtgtgtggctcggacgaagtgtttccagaccctgatctactgcatgtcctg
cctgtgtctgtgtcttgcagcagtcattgaccaaagctgtctacagcttgagagcctctgcaggccagggctgtctgtgcacatcctacgttactattcaaaactgcacagcagcatgaaaaacaggccctcttctctcatttacacatatgtgaagaagacacagcaagtgagaaaacgtg
accgaaagccacgcggacaggtggccgagggccaaacccaacttcctgtgatgtga