

RASL10B-RAS-like, family 10, member B Gene

Reference: PTXBC041133

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	alternative protein RASL10B; RRP17; VTS58635; ras-like protein family member 10B; Ras-related protein 17; ras-like protein VTS58635; RAS like family 10 member B		

For research use only.

GENE DESCRIPTION	RAS-like, family 10, member B
GENE ID	91608
ORIGIN SPECIES	Human
ACCESSIONS	BC041133
VECTOR	pDONR223
NCBI SYMBOL	RASL10B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtctccacctaccgggtgcccgtgctggggcgcgaggtgtgggcaagagtgccatcgtgcgccagttctgtacaagagttcagcgaggtctgctccccaccaccgccgcccgtttacctgctgctgctcatgaacggccacgtgcacgacctccagatcctgacttc
caccatcagcgcctccctgtcaatacgtccaggagtgggcgagacacctgctgcaggggactccggagtgccacgcctacatcctggtctacgacatctgctctttgacagctttgagtagctcaagaccatccgccagagatcctggagacgaggggtgatcggaacctcagaga
cgcccatcatcgtgggcaacaagcgggacctgcagcgcggacgctgatcccgctggaacgtgtgcacacctgttacgcaagacctggaagtgcggctacgtggaatgctcgccaagtacaactggcacatcctgctcttcagcagctgctcaagagcgtcggtgcg
cccgttgcaagcagctgcacgctgcccctgcgtccaggcgcgctgcgccgaaccgctgcgccatcatgtga