



Reference: PTXBC041374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NGL-1; NGL1; leucine-rich repeat-containing protein 4C; netrin-G1 ligand; leucine rich repeat containing 4C		

For research use only.

GENE DESCRIPTION	leucine rich repeat containing 4C
GENE ID	57689
ORIGIN SPECIES	Human
ACCESSIONS	BC041374
VECTOR	pDONR223
NCBI SYMBOL	LRRC4C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]