

RUNDC2B-RUN domain containing 2B Gene

Reference: PTXBC041583

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	RUN domain containing 2B
GENE ID	400509
ORIGIN SPECIES	Human
ACCESSIONS	BC041583
VECTOR	pDONR223
NCBI SYMBOL	RUNDC2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgaagagaggtccagatgcttctaccatggcagcaggtccgaactccatactctttgcgattaacattgacaacaaggattgaaacgggcagagtaagtgttgcctccaccgtttcagacctcttaaggaggtcaacgcagaatgtgacctgtgtaaggagtcacgcaaggagtg
agcagcgtgttcaggagatcacagcctcctctgccatctccatctcatcaaacctgaacaggagaccgacctgtccctgtgtgtccaggaatgtcagtgctgatgccaaatgcaaaaaggagcgggaagaagaaaaagcaagtgaccaacattatctcattgatgatgaggagatg
agcagaactctgggacatgtttaaaaagacacctggggcaggggagagctcagaggacaactccgaccactcctgttcaatatcatgtccgctttgaaagcccttcgggccaactccaatgggaagtgcagagcagcaactctgtgaaaattgattccctgtcttgaacaggaggttt
gggtaccagaagcttgatgtgaaaagcatcgatgatgaagatgtgatgaaaacgaagatgacgtgtatggaaactcatcaggacggaagcacagggccactcagagtcgcccagaaaaacggagcacactctgtcaccaggctggagtgcaatggcacgatctcggtcactg
caacctctgcctccgggtcaggtga