



Reference: PTXBC041611

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD158D; G9P; KIR-103AS; KIR-2DL4; KIR103; KIR103AS; killer cell immunoglobulin-like receptor 2DL4; CD158 antigen-like family member D; MHC class I NK cell receptor KIR103AS; killer cell immunoglobulin-like receptor, two domains, long cytoplasmic tail, 4;		

For research use only.

GENE DESCRIPTION	killer cell immunoglobulin-like receptor, two domains, long cytoplasmic tail, 4
GENE ID	3805
ORIGIN SPECIES	Human
ACCESSIONS	BC041611
VECTOR	pDONR223
NCBI SYMBOL	KIR2DL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]