

ZSCAN2-zinc finger and SCAN domain containing 2 Gene

Reference: PTXBC041620

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZFP29; ZNF854; zinc finger and SCAN domain-containing protein 2; zfp-29; zinc finger protein 29 homolog; zinc finger protein 854; zinc finger and SCAN domain containing 2		

For research use only.

GENE DESCRIPTION	zinc finger and SCAN domain containing 2
GENE ID	54993
ORIGIN SPECIES	Human
ACCESSIONS	BC041620
VECTOR	pDONR223
NCBI SYMBOL	ZSCAN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggctgcagacatcccagagtgaccactccgctgagctccttggccagggtgacctcaagaggaagatagacaggaggaggaggtcaccaccatgatcctggaggatgactcctgggtgcaagaagctgtgctgcaggaggatggccctgagctgagcccttcccagagtg
ctggcaaggcgccccagaggaggtgaccaggggaccacagggtgcactcgccgcctccgagagctctgccggcgctgctgagaccagaggtacaccaaggagcagatgttaacatgctgccaaaggaaatcaggcttggctgcaagagcatggcctgaaagca
gtgaggaggcagcgccctggtggaagactgaccagaccctcaggacagtgcgtgtagctgtcttcttcctccggtggaggtgaccagttgttaa