

## NPL-N-acetylneuraminate pyruvate lyase (dihydrodipicolinate synthase) Gene

Reference: PTXBC042003

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C112; C1orf13; NAL; NPL1; N-acetylneuraminate lyase; N-acetylneuraminate pyruvate lyase (dihydrodipicolinate synthase); N-acetylneuraminic acid aldolase; NALase; dihydrodipicolinate synthetase homolog 1; sialate-pyruvate lyase; sialic acid aldolase; N-ace |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | N-acetylneuraminate pyruvate lyase (dihydrodipicolinate synthase)                                 |
| <b>GENE ID</b>          | 80896                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC042003                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | NPL                                                                                               |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggccttccaaagaagaacttcagggtctgtgctgcaaccatcagccaatgactgagaatggagaatcaactttcagtaattggtcagtatgtggattatctgtgaaagaacaggaggagtgagaacatttttgaatggcacaacaggagaaggcctgtccctgagcgtctcag  
agcgtgccagggttcagaggagtggttgacaaaagggaaggacaagctggatcagggtgataatcacgtaggagcactgagcttgaaggagtcacaggaaactggccgaacatgcagcagaataaggagctgagggcatcgctgattgcaccgttcttctcaagccatggaccaa  
agatatctgattaattttctaaaggagtggtgctgcccctcgcctgccaattttattactatcacattcctgacagggtgaaagattcgtgctgaggagttgttgatgggattctggataagatccccacctccaagggtgaaatcagtgatacagatctcttagactcgggca  
atgtgtgatcagaatgccagcaacagtttcttcttttgggttgatgagcaattgttggtgctgctggtgaggaactggagcagtgggcagtcactataactacctgggaaaaagacaaccagatgttgaggcttttgaacaaaggactctcttag