

C3orf46-chromosome 3 open reading frame 46 Gene

Reference: PTXBC042038

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 3 open reading frame 46
GENE ID	255330
ORIGIN SPECIES	Human
ACCESSIONS	BC042038
VECTOR	pDONR223
NCBI SYMBOL	C3orf46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaagggcactgtgtctggcctgtgaggtgtggtgatgcagagactggcaaggtgtgcatcatctcaggacctcgttcaggtgaaggtgtgtgtctgttgagccgtgaggatccgagccccatcatgcagataaggatgggcaccagctgggaagaattggatacatctc
ccactgtgtgaaggaggcagcacagcttctggggcattccagtgctggcaccctgtgtcaccaggctgacttctccccacgatgtccatctatatcactggatcaccaaccaccagaacctttctctttggcaatgccatgccagcctgacctccactggtctgtcacaag
cgggacatcctggatctcaaggcagcatcacaagatgagccattaa