

MBNL3-muscleblind-like 3 (Drosophila) Gene

Reference: PTXBC042090

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHCR; MBLX; MBLX39; MBXL; muscleblind-like protein 3; Cys3His CCG1-required protein; muscleblind-like 3; muscleblind-like X-linked protein; muscleblind like splicing regulator 3		

For research use only.

GENE DESCRIPTION	muscleblind-like 3 (Drosophila)
GENE ID	55796
ORIGIN SPECIES	Human
ACCESSIONS	BC042090
VECTOR	pDONR223
NCBI SYMBOL	MBNL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgcgccagcagatgcagcttatgtccaaaacgtcaaatgtcatcacttggctctttctatgactccatcaattccagctaatcccatggcttcaatcctacataccacatcctggatgggctcgttctgcagaactgtaccaatacacctgttctgattcctggaaccacac
ctcttgcaatgccaggagctgttggcccaaaatgatgcttcagataaaactggaggtttgccgagaatttcagcgtggaattgtacccgtggggagaaatgattgccgctatgctcacccctactgatgcttccatgattgaagcagtgataactgtgacaatctcatggattacataaag
gtcgtatgctcggggagaaatgcaagtactttcatctcctgcacacttgcaagccagactcaaggcagctcatcatcagatgaaccattcagctgcctctgccatggccctgcagcctgggtacactgcaactgataccaaagagatcagcactggaaaagcccaatgggtgccacccgggt
ctttaatccactgttttccactgccaacaggctctgactaaactgcagctccacagccggcatttatccctgcagggccaatactgtgcatggcaccgcctcaaatattgtgccatgatgcacggtgctacacctaccactgtgtctgcagcaacaacacctgccaccagcgttccgttcg
ctgcaccaactacaggcaatcagctgaaattctga