



Reference: PTXBC042174

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CK-4; CK4; CYK4; WSN1; keratin, type II cytoskeletal 4; cytokeratin 4; keratin 4, type II; type-II keratin Kb4; keratin 4		

For research use only.

GENE DESCRIPTION	keratin 4
GENE ID	3851
ORIGIN SPECIES	Human
ACCESSIONS	BC042174
VECTOR	pDONR223
NCBI SYMBOL	KRT4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr