



Reference: PTXBC043165

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf146; FIGLER8; SALM5; leucine-rich repeat and fibronectin type-III domain-containing protein 5; fibronectin type III, immunoglobulin and leucine rich repeat domains 8; leucine rich repeat and fibronectin type III domain containing 5		

For research use only.

GENE DESCRIPTION	leucine rich repeat and fibronectin type III domain containing 5
GENE ID	145581
ORIGIN SPECIES	Human
ACCESSIONS	BC043165
VECTOR	pDONR223
NCBI SYMBOL	LRFN5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]