



Reference: PTXBC043277

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INPP5H; synaptojanin-2; inositol phosphate 5'-phosphatase 2; inositol polyphosphate-5-phosphatase H; phosphoinositide 5-phosphatase; synaptic inositol 1,4,5-trisphosphate 5-phosphatase 2; synaptojanin 2		

For research use only.

GENE DESCRIPTION	synaptojanin 2
GENE ID	8871
ORIGIN SPECIES	Human
ACCESSIONS	BC043277
VECTOR	pDONR223
NCBI SYMBOL	SYNJ2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]