



Reference: PTXBC043367

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM26A; bA225H22.7; calcium homeostasis modulator protein 3; family with sequence similarity 26, member A; calcium homeostasis modulator 3		

For research use only.

GENE DESCRIPTION	calcium homeostasis modulator 3
GENE ID	119395
ORIGIN SPECIES	Human
ACCESSIONS	BC043367
VECTOR	pDONR223
NCBI SYMBOL	CALHM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]