



Reference: PTXBC045535

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf25; MST070; MSTP070; TRM1L; bG120K12.3; TRMT1-like protein; TRM1 tRNA methyltransferase 1-like; TRM1-like protein; tRNA methyltransferase 1 homolog-like; tRNA methyltransferase 1 like		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 25
GENE ID	81627
ORIGIN SPECIES	Human
ACCESSIONS	BC045535
VECTOR	pDONR223
NCBI SYMBOL	C1orf25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]