

ZNF428-zinc finger protein 428 Gene

Reference: PTXBC045799

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf37; Zfp428; zinc finger protein 428; enzyme-like protein PIT13		
For research use only.			

GENE DESCRIPTION	zinc finger protein 428
GENE ID	126299
ORIGIN SPECIES	Human
ACCESSIONS	BC045799
VECTOR	pDONR223
NCBI SYMBOL	ZNF428
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagagacccgtgagccagctgagactgggggctacgccagcttggagaagatgatgatgaccttccccaggccccgagcattcctctgattcagaatacactctctcagagccggactccgaagagggaagaagatgaggaggagggaagaggagaccactgacatcct
gaatatgatcctggctacaaggtgaagcagcgccttggggggccgtgtgtgccccatccccggggccccccgtgcagccagccccgggccccagccttgcagctctgtggccctcacccttggggaggccccaccgggaacccccctgcccgtctgtgcccctgtac
agccccccaggaagcaccagccccctgaaggcagggccctcggggagggaagaggaggaaccacctcgggctggggaggggccgaccagctgggccccggaggaggagggaagaggaggagggaacctaccactgtacggaatgtgaggattcctcgacaacctggggg
agctgcacgggcacttcatgctgcatgccccgggtgaggtgtag