

DNALI1-dynein, axonemal, light intermediate chain 1 Gene

Reference: PTXBC046117

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P28; dJ423B22.5; hp28; axonemal dynein light intermediate polypeptide 1; dJ423B22.5 (axonemal dynein light chain (hp28)); dynein, axonemal, light intermediate polypeptide 1; inner dynein arm light chain, axonemal; inner dynein arm, homolog of clamydomonas		

For research use only.

GENE DESCRIPTION	dynein, axonemal, light intermediate chain 1
GENE ID	7802
ORIGIN SPECIES	Human
ACCESSIONS	BC046117
VECTOR	pDONR223
NCBI SYMBOL	DNALI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattccgcccgcagactcttgcctcaagtagacacccagtgctggtgagccggaacacggagaacacggagcccaaggtcggctactgaaagtcagccccagcagcctggacctcaggttcagccccacagccaccaagaccaagctcccctcaactccctgtgtccag
atcctacaaagcaggcagaagaatcttgaatgccatactaccccaaggagtggtggaagacacgcagctatggatccagcaggtgtccagcacccctagcaccaggtggagctggtgacctccaggagcagttagacttaagctgcagcagcggcaggccagggaaca
ggcatctgcccgtccgagggaactctactcacagtgtttgatgagttgatccggaggtcacatcaactgtgcgagagggggctgctgctgctgcgaggtccgggagagatccgcatgacatcgctgcctaccagacctgtacgagagcagcgtggcgtttggcatgaggaag
gcactgcaggctgagcaggggaagtcagacatggagaggaaatcgagaattggagacggaaaagagagacctggagagcgaagtgaacgagcagaaggcaaatgtgaagccactgagaagcgggagagcagagggcgcaggtggaggagaagaagcacaatgagga
gattcagttcctgaagcgaacaatcagcagctgaaggcccaactggaagcattattgcaccaagaagtga