

THAP8-THAP domain containing 8 Gene

Reference: PTXBC046207

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	THAP domain-containing protein 8; THAP domain containing 8		
For research use only.			

GENE DESCRIPTION	THAP domain containing 8
GENE ID	199745
ORIGIN SPECIES	Human
ACCESSIONS	BC046207
VECTOR	pDONR223
NCBI SYMBOL	THAP8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagtactcaggcgccgaactgctccaactgcggcgccctgggtgcagacaaccgcccgtgagcttctacaagttccactgaaggatgtcccccggctgcaggcctggctgcagcacatgggctgtgagcactgggtgccagctgccaccagcactgtgcagcg
agcacttcacaccctcctgtctccagtggcgctgggtgtgctactgcggcctgatgcagtgcctccatcttctccggggaccacctgccaaagatcagcggaggacccgaagcaccagaagccagtctgcgcgcgcctccctacagaagaatcacccctgccccagag
cctgcatcccatctctgtggccagtgcgcctagtgtgtgtggccccacatcggggagcggcgaagactgtgcccacatgctcctgacccccctggccccctgcgccaacttctgagcgggtcacaacctgaagtcctgcccacagggccagaccggctggccagctgtggg
agcactgcaacgcgggtgcggaggctgcaacgggtgccaggagcggcaccagggcgagctgcaggccctggaacgggtggcagcagctacacggggagagcctgtgtggcacgggcacgcgggtctgcagcgctgacaacagcccagacccttgacctgaggaatcc
caaaccttcaccatcatctgtggaggcctgacatagccatgtcttggccaggaccctgcacctgccacagtgatgccaagccggagctcctggacactcggatccccagtgataa