

STX12-syntaxin 12 Gene

Reference: PTXBC046999

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	STX13; STX14; syntaxin-12; syntaxin 12		
For research use only.			

GENE DESCRIPTION	syntaxin 12
GENE ID	23673
ORIGIN SPECIES	Human
ACCESSIONS	BC046999
VECTOR	pDONR223
NCBI SYMBOL	STX12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcatcgggtcccttagacatgtaccggaacccggggccctcggggcccccagctccgggaactcagcagcatcatccagacgtgcagcggcaacatccagcggatcagccaagccactgctcagataaagaatttgatgagccagctaggaactaagcaggactcaagcaagctac
aggaaaaatctgcaacagttacacactccacaatcagctcgccaaggaaacaaatgaattgctgaaagaatttagggctctgcccttcccttatctactcagaacagcggccagcagagacttcagaaggaacgctcatgaatgacttctctgagccttaacaatttcaggctgtgca
gagaagggtatctgaaaaggaaaaggagagattgcccagagcaagagctggaatctctcttctgagaagagagggcaagagagagcagctgggtctcatgtgacagccatgaggagtggaaccagatgcagagccaggagaggtggccatcactgagcaggatttggaaactt
attaaagaaagagaaacggcaattcggcagctggaggctgacatttggatgtcaatcagatatttaaagatttggccatgatgatccatgaccagggtgatctgattgatagcatagaagccaatgtggaaagctcagaggtgcacgtcgaagagccactgaacagttacagcagctgc
ttactatcagaaaaatctcgcaagaagatgtgtatcctggtgcttgcctgtcagtgattattctaattctgggaacttattatctgctagtttataaaacgaagtga