



Reference: PTXBC047221

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SEC63 homolog, protein translocation regulator; SEC63-like protein; SEC63-like; SEC63 protein translocation regulator; SEC63 homolog; translocation protein SEC63 homolog; DNAJC23; ERdj2; PCLD2; PRO2507; SEC63L		

For research use only.

GENE DESCRIPTION	SEC63 homolog (<i>S. cerevisiae</i>)
GENE ID	11231
ORIGIN SPECIES	Human
ACCESSIONS	BC047221
VECTOR	pDONR223
NCBI SYMBOL	SEC63
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]